




BOSS
INSTRUMENTS LTD





IMAGE-TRAC™
POSTERIOR SPINE RETRACTOR SYSTEM

RETRACTOR FRAMES


The **BOSS Image-Trac™ Retractor System** provides superior exposure for all surgical procedures from T1 to S1.

Lumbar Fusion (single or multi-level)
Posterior Lumbar Interbody Fusion (PLIF)
Hemi-Laminectomy
Lumbar Laminectomy
Micro-Lumbar
Discectomy
Spinal Fusion
for Scoliosis



Scoville Frame
73-8010

Lumbar Frame
73-8012



RETRACTOR
TRAY



RETRACTOR BLADES

The **BOSS Image-Trac™ Retractor System** features simplicity as well as versatility in posterior lumbar retraction.


Titanium Retractor Blades:

Assorted widths and lengths to meet the needs of all patient anatomies
Designed to allow for maximum visualization
Radiolucent to permit intraoperative imaging



Quick-attach blade holder #73-8018 provides unsurpassed flexibility in blade attachment anywhere on the lumbar retractor frame

Micro-adjust blade holder #73-8016 permits delicate retraction in the difficult L4 and L5 areas



**NARROW
BLADE
TRAY**



**WIDE
BLADE
TRAY**

ORDERING INFORMATION

Image-Trac™ Complete Set - 73-8000

Item	Description	Units in set
ACCESSORIES		
73-6095	Retractor Handle, Side-Loading	2
73-8001	Case for Complete Image-Trac	1
73-8010	Scoville Body, Double-hinged, 200mm Rack	1
73-8012	Lumbar Retractor Body for Multiple Blades	1
73-8014	Lumbar Retractor T-Handle for 73-8012	1
73-8016	Lumbar Retractor Adjustable Blade Holder for 73-8012	4
73-8018	Lumbar Retractor Blade Holder for 73-8012	4
TITANIUM RETRACTOR BLADES		
73-8040T	Sacral Cutout Blade, Right, 15 deg, w/teeth, 30x40mm	1
73-8042T	Sacral Cutout Blade, Left, 15 deg, w/teeth, 30x40mm	1
73-8044T	Sacral Cutout Blade, Right, 15 deg, w/teeth, 30x45mm	1
73-8046T	Sacral Cutout Blade, Left, 15 deg, w/teeth, 30x45mm	1
73-8051T	Angled Out Blade 15 deg, w/teeth 45x50mm	2
73-8052T	Angled Out Blade 15 deg, w/teeth 45x60mm	2
73-8053T	Angled Out Blade 15 deg, w/teeth 45x70mm	2
73-8054T	Angled Out Blade 15 deg, w/teeth 45x80mm	2
73-8055T	Angled Out Blade 15 deg, w/teeth 60x50mm	2
73-8056T	Angled Out Blade 15 deg, w/teeth 60x60mm	2
73-8057T	Angled Out Blade 15 deg, w/teeth 60x70mm	2
73-8058T	Angled Out Blade 15 deg, w/teeth 60x80mm	2
73-8059T	Angled Out Blade 15 deg, w/teeth 60x100mm	2
73-8060T	Angled Out Blade 15 deg, w/teeth 60x120mm	2
73-8071T	Angled Out Blade 15 deg, w/teeth 23x35mm	4
73-8072T	Angled Out Blade 15 deg, w/teeth 23x45mm	4
73-8073T	Angled Out Blade 15 deg, w/teeth 23x55mm	4
73-8074T	Angled Out Blade 15 deg, w/teeth 23x65mm	4
73-8075T	Angled Out Blade 15 deg, w/teeth 23x75mm	4
73-8076T	Angled Out Blade 15 deg, w/teeth 23x85mm	4
73-8077T	Angled Out Blade 15 deg, w/teeth 23x100mm	4
73-8078T	Angled Out Blade 15 deg, w/teeth 23x120mm	4
STAINLESS STEEL RETRACTORS		
73-8022	Hook, shaft w/teeth sharp, 40mm	1
73-8023	Hook, shaft w/teeth sharp, 50mm	1
73-8024	Hook, shaft w/teeth sharp, 60mm	1
73-8025	Hook, shaft w/teeth sharp, 70mm	1
73-8030	Lumbar Spine T-Bar, 25x40mm	1
73-8032	Lumbar Spine T-Bar, 25x50mm	1
73-8034	Lumbar Spine T-Bar, 25x60mm	1
73-8036	Lumbar Spine T-Bar, 25x70mm	1

Item	Description
OPTIONAL RETRACTORS	
73-8026	Single Hook, Sharp, 40mm
73-8027	Single Hook, Sharp, 50mm
73-8028	Single Hook, Sharp, 60mm
73-8029	Single Hook, Sharp, 70mm
73-8029X	Single Hook, Sharp, 80mm
73-8063T	Angled Out Blade 15 deg, w/teeth 60x140mm



Image-Trac™ complete in a sealed container - 73-8008



Image-Trac™ complete in a sterilization case - 73-8000

BOSS
INSTRUMENTS LTD

104 Sommerfield Drive
Gordonsville, VA 22942 USA
www.bossinst.com

Toll Free: 800-210-2677
Phone: 540-832-5000

Email: inquiry@bossinst.com
Fax: 540-832-5515