




BOSS

INSTRUMENTS LTD

**SINUS
INSTRUMENTS**




"CHAIN™" SINUS FORCEPS		
length	bite	item
4½" (11.5cm)	2.5mm wide x 2.0mm high	93-6080
	3.0mm wide x 3.5mm high	93-6082
	5.5mm wide x 5.0mm high	93-6084*



CARLEY ROTATING BACKBITER		
length	curved	item
3½" (9.5cm)	2mm x 4mm, down	93-5698*
	2mm x 4mm, up	93-5699

MICRO ROTATING BACKBITER		
length	description	item
4½" (11.5cm)	5mm x 1.5mm bite	93-5634

THRUBITE™ FORCEPS		
length	description	item
5" (13cm)	3mm x 6mm bite, horizontal jaws	93-5730*
	3mm x 6mm bite, vertical jaws	93-5732






TOBEY ROTATING MICRO BACKBITER

length	description	item
4½" (11.5cm)	small oval shaped 3mm x 2mm jaws, opens greater than 90°	93-5638



THRUBITE™ PUNCH FORCEPS

length	description	item
3½" (9.5cm)	pediatric 1.5mm x 6.0mm bite, straight	93-7082
	pediatric 1.5mm x 6.0mm bite, 45°	93-7084*
	size 00, 2.5mm x 7.0mm bite, 45°	93-7090
	size 00, 2.5mm x 7.0mm bite, straight	93-7092
	size 0, 3.0mm x 8.0mm bite, straight	93-7094
	size 0, 3.0mm x 8.0mm bite, 45°	93-7096



"THRUBITE™" FRONTAL SINUS GIRAFFE FORCEPS

length	description	item
6" (15cm)	45°, 2mm vertical jaws	93-7050
	45°, 2mm horizontal jaws	93-7051
	45°, 3mm vertical jaws	93-7054
	45°, 3mm horizontal jaws	93-7055*
7" (18cm)	90°, 2mm vertical jaws	93-7052
	90°, 2mm horizontal jaws	93-7053



MURPHY "SLIM" BACKBITER		
length	description	item
3½" (9.5cm)	small oval shaped 3mm x 2mm bite, curved up	93-5636*
	small oval shaped 3mm x 2mm bite, curved down	93-5637



CIRCULAR CUTTING PUNCH		
length	description	item
7" (18cm)	straight, 3.5mm	93-6060*
	up 45°, 3.5mm	93-6061
	straight, 4.5mm	93-6065
	up 45°, 4.5mm	93-6066

SILVER SEPTAL NEEDLE HOLDER	
length	item
2½" (6.5cm)	93-7001
4½" (11.5cm)	93-7000*






SIDEBITING ANTRUM PUNCH

length	description	item
4" (10cm)	sidebiting antrum punch, right	93-5620*
	sidebiting antrum punch, left	93-5621

ROTATING ANTRUM PUNCH

length	description	item
4" (10cm)	BOSS rotatable antrum punch, pediatric	93-5640*
4½" (11.5cm)	BOSS rotatable antrum punch, adult	93-5641



"THRUBITE™" WEIL-BLAKESLEY FORCEPS

length	description	item
3½" (9.5cm)	straight, size 00, 2.5mm	93-6000
	straight, size 0, 3.5mm	93-6001
	straight, size 1, 4.0mm	93-6002
	45°, size 00, 2.5mm	93-6010
	45°, size 0, 3.5mm	93-6011
	45°, size 1, 4.0mm	93-6012
	angled right, size 0, 3.5mm	93-6015
	angled left, size 0, 3.5mm	93-6016

* indicates product pictured



INNER-FLOW™ GRUENWALD PUNCH FORCEPS

length	description	item
4" (10cm)	straight, 2mm jaws	93-5790*
	up, 2mm jaws	93-5792

SOFT TISSUE PEAR FORCEPS


length	description	item
4" (10cm)	upturned 45°	93-5786*
	upturned 80°	93-5787

INNER-FLOW™ NASAL SUCTION FORCEPS

length	description	item
4" (10cm)	straight shaft, straight jaws, 2.5mm	93-5920
	straight shaft, straight jaws, 3.5mm	93-5940*
	straight shaft, straight jaws, 4.0mm	93-5950
	straight shaft, angled jaws, 2.5mm	93-5960
	straight shaft, angled jaws, 3.5mm	93-5980
	straight shaft, angled jaws, 4.0mm	93-5990



KUHN FRONTAL SINUS SURGERY SET		
length	description	item
6" (15cm)	giraffe forceps, 45°, 2mm vertical jaw	93-5750
	giraffe forceps, 45°, 2mm horizontal jaw	93-5753
	giraffe forceps, 90°, 2mm horizontal jaw	93-5756
	giraffe forceps, 90°, 2mm vertical jaw	93-5759
	giraffe forceps, 45°, 3mm vertical jaw	93-5762
	giraffe forceps, 45°, 3mm horizontal jaw	93-5765
7" (18cm)	curette, small oblong, forward cutting, 90°	93-5548
7½" (19.5cm)	curette, small oblong, forward cutting, 45°	93-5545
	ostium seeker double-ended	93-5552



* indicates product pictured



FRONTAL SINUS INSTRUMENTS		
length	description	item
7" (18cm)	ostium probe	93-5770
	prograde rasp, coarse tip, 3mm	93-5775
	prograde rasp, fine tip, 3mm	93-5776
	curette, curved 75°	93-5778
	curette, curved 100°	93-5779
7½" (19.5cm)	suction tube, curved, 3mm	93-5772
	prograde rasp with suction channel, 4mm	93-5774

BOSS
INSTRUMENTS LTD

104 Sommerfield Drive
Gordonsville, VA 22942 USA

Toll Free: 800-210-2677

Phone: 540-832-5000

Fax: 540-832-5515

Email: inquiry@bossinst.com

www.bossinst.com