

ENT Theatres & Clinics



Precision Surgical has a growing range of ENT products for both Theatres and Clinics

- LED Headlights
- ENT Instruments
- Hand-held LED Lightsource
- Single Use Suction Tubes

Quantum Nano Headlight under 150g (including Li-ion battery)



Headband

- Wireless and compact, ingenious and futuristic design
- Innovative polymer material - very durable and does not deform
- Ergonomic design ensuring comfortable fit

Features

- Battery Indicator - 5 levels (20% increments)
- Light Intensity Control - 5 levels
- Supplied with 2 x Li-ion Batteries, Charger & Transport Case

Cable Hose & Housing

- Focussed beam can be adjusted 360°
- Hose retains its position even after years of use

Light Performance

- Shadow free light from the smallest medical housing
- Intelligent constant-current control to achieve stable illumination



Clinical Headlight - Black (PS3000-00)

Forward-facing LED Housing 40,000 Lux at 200mm

Surgical Headlight - Blue (PS4000-00)

Downward-facing LED Housing 50,000 Lux at 200mm

Designed by the world's leading optical experts and manufactured by German nano-scale processes, with the precise light-averaging system, the output light spot is powerful and perfectly defined. The sophisticated design significantly reduces the housing size and realises excellent coaxial and shadow-free lighting with the smallest housing diameter in the industry.

ENT Clinic Instruments

High quality range of reusable ENT instruments, commonly used in most ENT Outpatients Departments


Speculae

- Tumarkin
- Gruber
- Hartmann
- Thudicum



Micro Scissors & Forceps

- Bellucci Style Scissors
- Micro Cup Forceps
- Crocodile Style Forceps
small / medium / large
- Tilley Nasal Forceps




Probes & Hooks

- Jobson Horne
- St Barts Blunt Hook
- Formby Scoop & Hook
- Cawthorne Hooks



Other Clinic Instruments

- Zoellner Suction Tubes & fine ends
- Lacks Tongue Depressor
- Iris Scissors, straight & curved



Other Clinic Instruments available please contact for further details

ENT Theatre Instruments

An extensive high quality range of speciality ENT instruments, commonly used in all ENT Theatres

Micro Head & Neck

- Acland Clamps
- Micro Surgical Instruments




General & Micro Ear

- Bellucci Style Scissors
- Micro Cup Forceps
- Crocodile Style Forceps
small / medium / large



Nasal & Rhinoplasty

- Specialist Scissors
- Bottle & Killian Speculae
- Rasps, Elevators, Chisels



Sinuscopy

- Specialist Punches
- Sinus Forceps, Scissors
- Sinusscopes



Micro Laryngoscopy

- Laryngoscopes
- Micro Laryngeal Instruments
Kleinsasser



Extensive range of Theatre Instruments available please contact for further details

LED Portable Lightsource

Rechargeable battery powered LED light source which can be attached to any ENT scope or light guide carrier. No mains powered light source or fibre light cables are required. Can also be used for Laryngoscopy in Theatres



for Clinic and
Theatre use

Li-ion Clinic LED Portable
Lightsource (PS5000-00)

including 2 x Li-ion batteries, charger & pouch



light guide carrier
for Laryngoscopy

NiMH Theatre LED Portable
Lightsource (PS6000-00)

including 2 x NiMH batteries, charger & pouch

- Natural Colour Temperature 5,000k
- Light Intensity 1 - 570 Lumens
- Weight 80g
- Dimensions 100mm length x 23mm diameter
- Operating time - 45mins at 500+ lumens (Li-ion)
- Operating time - 120mins at 225+ lumens (NiMH)
- Recharge time - 2-3 hours

Rigid Telescopes (manufactured in Germany) available
please contact for further details

