


Single Use Suction Tubes for ENT Surgery



- ▶ Strong Surgical Steel cannula for dissection
- ▶ Machine manufactured & 100% quality guarantee
- ▶ Rounded atraumatic distal tips
- ▶ Internal uniform diameter creates a "non blocking design"
- ▶ Single sterile blister packaging (gamma irradiation)
- ▶ Colour coded & well-balanced handles



Laryngoscopy

Code	Internal Diameter mm	External Diameter mm	fg	*Cannula Length mm	Colour	Box
CSSR20X330	1.45	2mm	2fg	330	Green	20
CSSR30X300	2.54	3mm	3fg	330	Blue	20

Nasal

Adson Frazier	Adson Frazier with Stilette	Internal Diameter mm	External Diameter mm	Working Length mm	fg	Colour	Box
CSSR20X80	CSSR20X80S	1.45	2mm	70mm	6fg	Green	60
CSSR30X80	CSSR30X80S	2.54	3mm	70mm	9fg	Blue	60
CSSR40X80	CSSR40X80S	3.0	4mm	70mm	12fg	White	60
CSSR20X130	CSSR20X130S	1.45	2mm	100mm	6fg	Green	60
CSSR30X130	CSSR30X130S	2.54	3mm	100mm	9fg	Blue	60
CSSR40X130	CSSR40X130S	3.0	4mm	100mm	12fg	White	60
CSSR20X160	CSSR20X160S	1.45	2mm	130mm	6fg	Green	60
CSSR30X160	CSSR30X160S	2.54	3mm	130mm	9fg	Blue	60
CSSR40X160	CSSR40X160S	3.0	4mm	130mm	12fg	White	60

Ear

Verhoeven House Suction Tubes	Internal Diameter mm	External Diameter mm	swg	fg	*Cannula Length mm	Working Length mm	Colour	Box
CSSR07X70	0.43	0.7	24	2fg	70	60	Black	60
CSSR10X70	0.7	1	22	3fg	70	60	Brown	60
CSSR14X70	0.96	1.4	20	4fg	70	60	White	60
CSSR07X90	0.43	0.7	24	2fg	90	80	Black	60
CSSR10X90	0.7	1	22	3fg	90	80	Brown	60
CSSR14X90	0.96	1.4	20	4fg	90	80	White	60