


# Single Use Suction Tubes for Spinal & Neurosurgery



- ▶ Strong Surgical Steel cannula for dissection
- ▶ Machine manufactured & 100% quality guarantee
- ▶ Rounded atraumatic distal tips
- ▶ Internal uniform diameter creates a “non blocking design
- ▶ Single sterile blister packaging (gamma irradiation)
- ▶ Colour coded & well-balanced handles
- ▶ Available with stilette



Adson Frazier	Adson Frazier with Stilette	Internal Diameter mm	External Diameter mm	Working Length mm	fg	Colour	Box
CSSR20X80	CSSR20X80S	1.45	2mm	70mm	6fg	Green	60
CSSR30X80	CSSR30X80S	2.54	3mm	70mm	9fg	Blue	60
CSSR40X80	CSSR40X80S	3.0	4mm	70mm	12fg	White	60
CSSR20X130	CSSR20X130S	1.45	2mm	100mm	6fg	Green	60
CSSR30X130	CSSR30X130S	2.54	3mm	100mm	9fg	Blue	60
CSSR40X130	CSSR40X130S	3.0	4mm	100mm	12fg	White	60
CSSR20X160	CSSR20X160S	1.45	2mm	130mm	6fg	Green	60
CSSR30X160	CSSR30X160S	2.54	3mm	130mm	9fg	Blue	60
CSSR40X160	CSSR40X160S	3.0	4mm	130mm	12fg	White	60