

Trimline Style ACDF Set

Aluminium, Colour Coded Blades

PRECISION SURGICAL
Instruments & Retractors of Distinction



**Stealth Side Loading
Transverse Retractor,
double hinge. 150mm rack**

SS1101



**Stealth Side Loading Transverse
Offset (step) Retractor, double
hinge, 150mm rack**

SS1102



Blunt Colour Blades (20mm), aluminium

SS1010	30mm depth	Silver
SS1012	40mm depth	Red
SS1014	50mm depth	Violet
SS1016	60mm depth	Gold
SS1018	70mm depth	Green




Toothed Medial Blades (23mm), aluminium

SS3531	30mm depth	Silver
SS3533	40mm depth	Red
SS3535	50mm depth	Violet
SS3537	60mm depth	Gold
SS3539	70mm depth	Green



Trimline Style ACDF Set


Aluminium, Colour Coded Blades



**Stealth Blade
Loading Handle**

SS0806

Qty	Code	Recommended Trim-Line ACF Set
1	73-6104	Sterilisation Trays For Blades & Retractors
2	SS0806	Stealth Blade Loading Locking Retractor Handle
1	SS1101	Stealth Side Loading Transverse Retractor Double Hinge
1	SS1102	Transverse Offset Retractor
		20mm Wide Blunt Colour Coded Blades
2	SS1010	Stealth ACF Blade Blunt 30mm Silver
2	SS1012	Stealth ACF Blade Blunt 40mm Red
2	SS1014	Stealth ACF Blade Blunt 50mm Violet
2	SS1016	Stealth ACF Blade Blunt 60mm Gold
2	SS1018	Stealth ACF Blade Blunt 70mm Green
		23mm Wide Medial Colour Coded Blades
2	SS3531	Stealth Notched Medial Retractor Blade 23x30mm Silver
2	SS3533	Stealth Notched Medial Retractor Blade 23x40mm Red
2	SS3535T	Stealth Notched Medial Retractor Blade 23x50mm Violet
2	SS3537T	Stealth Notched Medial Retractor Blade 23x60mm Gold
2	SS3539T	Stealth Notched Medial Retractor Blade 23x70mm Green



**Sterilisation Tray for blades
& Retractor Bodies (with
silicone mat)**

73-6104

To customise your
Trimline set, please
contact us directly