



Guide your hand with cool, bright light.

LIGHTMAT® SURGICAL ILLUMINATION

**Deep cavity illumination
using your existing retractors**

BOSS
INSTRUMENTS LTD

104 Sommerfield Drive
Gordonsville, VA 22942 USA

Toll Free: 800-210-2677

Phone: 540-832-5000

Fax: 540-832-5515

Email: inquiry@bossinst.com

www.bossinst.com

LightMat® delivers cool, bright light deep into the surgical cavity, by attaching to existing retractors or other flat-surfaced devices. The patented fiber-optic technology provides white, shadowless light throughout the procedure, often eliminating the need for headlights or costly new retractors. LightMat® is a cost-effective option for providing direct light, where you need it, when you need it.

UA2550 | TA2334

MM1462

OM2130

OM1242

LightMat


Sold 5 per box in individual, sterile packages.

UA2550	Ultra-Thin LightMat w/adhesive backing
TA2334	Ultra-Thin LightMat w/built-in malleable retractor
MM1462*	MicroBright LightMat w built-in malleable retractor
OM2130*	OrthoLight Spotlight w/adhesive backing
OM1242*	OrthoLight Spotlight w/ultra-malleable backing and optional adhesive strip

* Requires a 2.5 mm cable.

Reusable Fiber-Optic Cables

	4mm cable	2.5mm (8ft) cable	2.5mm (12ft) long cable
ACMI	#005011	#006241	#007986
Olympus	#005365	#006424	#008219
Storz	#005366	#006422	#008217
Wolf	#005757	#006423	#008218



Light Source

008554	50W BRICK Light Source w/power supply (US only)
--------	---

BOSS
INSTRUMENTS LTD