



Enlightened versatility.

VERSALIGHT™ SURGICAL ILLUMINATION

**Multiple Functions
and Illumination in a
Single-use, Handheld Device**

BOSS
INSTRUMENTS LTD

104 Sommerfield Drive
Gordonsville, VA 22942 USA

Toll Free: 800-210-2677

Phone: 540-832-5000

Fax: 540-832-5515


Email: inquiry@bossinst.com

www.bossinst.com

VersaLight is a multipurpose tool that can:

- Illuminate
- Irrigate
- Suction
- Provide Moderate Retraction

Available with or without a Frazier tip, the VersaLight uses patented MicroLens technology to direct cool, white light directly where you need it.



VersaLight

Sold 5 per box in individual, sterile packages.
Includes integrated suction/irrigation tubing.

VSG411	Mini VersaLight
VSX433	Mini VersaLight with extended Frazier tip
VSG311	VersaLight
VSX333	VersaLight with extended Frazier tip
VSG211	VersaLight without integrated tubing (not pictured)

Reusable Fiber-Optic Cables

	2.5mm (12ft) long cable
ACMI	#007986
Olympus	#008219
Storz	#008217
Wolf	#008218



Light Source

008554	50W BRICK Light Source w/power supply (US only)
--------	--

BOSS
INSTRUMENTS LTD

© 2018 BOSS Instruments, LTD. REV B L-65 BOSS and
BOSS Instruments, LTD. are registered trademarks.