

CASPAR CERVICAL

SPINE DISTRACTOR SETS



Vertebral Body Distractor
With Regular Length Toothed Bar

73-2588*
right

73-2587
left



Vertebral Body Distractor
With Elongated Toothed Bar
For Multilevel Distraction

73-2598
right

73-2597*
left



Distractor Drill Guides
For Parallel Positioning
of the Distraction Pins

73-2586*
right
used with:
73-2588 or 73-2587

73-2585
left
used with:
73-2587 or 73-2588



Distractor Twist Drill
Tip Diameter 1.7mm

73-2596*
8mm ($\frac{5}{16}$ ") depth
5 $\frac{3}{4}$ " (14.6cm)
used for pre-boring
holes for distraction



Distractor Screwdriver

73-2590*
8" (20cm)

