


Chisel Set of 6



Material Used: Aluminum Handle, Stainless Steel (AISI-420)

Total Length: 10³/₄" (270mm),

Handle Length: 5³/₄"(145mm), **Handle Diameter:** 28 mm, **Working Length:** 5" (125mm)

Available as a set, or individually on demand.



Straight

Curved

	Chisel	Sizes
Length	Overall: 10 ³ / ₄ " – 270mm Handle: 5 ³ / ₄ " – 130mm Working Point: 5" – 150mm	
Straight	PS-C-106	6mm
Straight	PS-C-110	10mm
Straight	PS-C-115	15mm
Straight	PS-C-120	20mm
Straight	PS-C-125	25mm
Straight	PS-C-130	30mm
Curved	PS-C-206	6mm
Curved	PS-C-210	10mm
Curved	PS-C-215	15mm
Curved	PS-C-220	20mm
Curved	PS-C-225	25mm
Curved	PS-C-230	30mm