

## Cobb Elevators



**Material Used:** Titanium / Stainless Steel (AISI-420)

**Total Length:** 12<sup>3</sup>/<sub>5</sub>" (320mm),

**Handle Length:** 7<sup>4</sup>/<sub>5</sub>" (198 mm), **Handle Diameter:** 28 mm, **Working Length:** 4<sup>7</sup>/<sub>8</sub>" (120mm)

Available as set, or individually in Titanium or Stainless Steel

Code	Description	Sizes	Quantity
PS-CE-1410	Cobb Elevator	10mm	1
PS-CE-1415	Cobb Elevator	15mm	1
PS-CE-1420	Cobb Elevator	20mm	1
PS-CE-1425	Cobb Elevator	25mm	1
PS-CE-1430	Cobb Elevator	30mm	1
PS-CE-1435	Cobb Elevator	35mm	1
PS-CE-1400	Elevator Rack		1