

Hardy Type Bayonet Instruments Set



Material Used: Titanium / Stainless Steel (AISI-420)

Titanium Coating: Blue / Black / Gold Total

Length: 9½", 240mm

Available as set, or individually in titanium or stainless steel

Code	Description	Size	Quantity
PS-HCU-1102	Hardy Type Curette, Bayonet, Angled Up	2mm	1
PS-HCU-1103	Hardy Type Curette, Bayonet, Angled Up	3mm	1
PS-HCU-1105	Hardy Type Curette, Bayonet, Angled Up	5mm	1
PS-HCD-1102	Hardy Type Curette, Bayonet, Angled Down	2mm	1
PS-HCD-1103	Hardy Type Curette, Bayonet, Angled Down	3mm	1
PS-HCD-1105	Hardy Type Curette, Bayonet, Angled Down	5mm	1
PS-HCL-1102	Hardy Type Curette, Bayonet, Angled 90° Left	2mm	1
PS-HCL-1103	Hardy Type Curette, Bayonet, Angled 90° Left	3mm	1
PS-HCL-1105	Hardy Type Curette, Bayonet, Angled 90° Left	5mm	1
PS-HCR-1102	Hardy Type Curette, Bayonet, Angled 45° Right	2mm	1
PS-HCR-1103	Hardy Type Curette, Bayonet, Angled 45° Right	3mm	1
PS-HCR-1105	Hardy Type Curette, Bayonet, Angled 45° Right	5mm	1
PS-HDR-1105	Hardy Type Curette, Bayonet, Distal Tip Curved Right	5mm	1
PS-HDL-1105	Hardy Type Curette, Bayonet, Distal Tip Curved Left	5mm	1

Hardy Type Bayonet Instruments Set

Code	Description	Size	Quantity
PS-HDU-1121	Hardy Type Dissector, Bayonet, Angled Up		1
PS-HDD-1131	Hardy Type Dissector, Bayonet, Angled Down		1
PS-HEL-1141	Hardy Type Enucleator, Bayonet, Micro Left		1
PS-HER-1142	Hardy Type Enucleator, Bayonet, Micro Right		1
PS-HEL-1145	Hardy Type Enucleator, Bayonet, Standard Left		1
PS-HER-1146	Hardy Type Enucleator, Bayonet, Standard Right		1
PS-HKV-1151	Hardy Type Knife Handle, Bayonet, Vertical (Not Shown)		1
PS-HKL-1152	Hardy Type Knife Handle, Bayonet, Lateral (Not Shown)		1
PS-HF-1161	Hardy Type Fork, Bayonet, 3 Prong		1
PS-HPS-1174	Hardy Type Pituitary Spoon, Bayonet, Curved Up	4mm	1
PS-HT-1100	Sterilizing Case for Hardy Type Bayonet Curettes		1