

Hintermann Distractors

Our range of specialist Distractors for Hand, Foot & Ankle Surgery



Hintermann Distractors with Open and Closed Arms

Hintermann Distractors are used in Foot and Ankle Surgery, providing flexible positioning without compromising the procedure.

Weinraub Spreader

The **Weinraub Spreader** helps prevent slippage or falling out of the joint by placing the arms on either side of the area to be distracted, driving two pins and opening the joint. The K-Wire size is 1.6mm.

We supply these available with 1.6mm and 2.5mm K-Wire holes, with open or closed arms, in either 15cm or 21cm lengths.

Hintermann Distractor with Closed Arms

Length 150mm
Length 210mm



Product Code	Length (mm)
PSB52-840	150
PSB52-845	210

Hintermann Distractor with Open Arms


Length 150mm
Length 210mm



Product Code	Length (mm)
PSB52-860	150
PSB52-865	210

Weinraub Joint and Calcaneal Spreader

Length 210mm



Product Code	Length (mm)
PSB52-900	210