



**BOSS**  
INSTRUMENTS LTD

**KERRISON & IVD  
RONGEURS**


# The Anatomy of a BOSS Kerrison Rongeur

## Spurling Style

The original handle style for all kerrison rongeurs, Spurling kerrison rongeurs feature a 15mm bite opening and a long handle reach.

Approximate dimensions:

Jaw Opening (A) - 15 millimeters  
Handle Width (B) - 140 millimeters  
Handle Height (C) - 145 millimeters



## Performance Coating


Boss offers all kerrison rongeurs in high-quality hardened German stainless steel. As kerrison rongeurs are used so frequently, and high performance is needed at all times, Boss has 2 performance coatings available to make your kerrison rongeurs keep their edge and smooth action for a longer period of time.

**P** = Panther™

A black ceramic coating is bonded to the steel at high heat and pressure, Panther rongeurs feature a matte black finish, smooth cutting action, and increased cutting edge hardness.

**SC** = Super-Coat

Super-Coat utilizes the strength and lubricity of Titanium Nitride bonded to the instrument surface. This special coating gives the rongeur its distinctive gold color.




### Micro Style

The most popular style of kerrison rongeur for both lumbar and cervical procedures. Micro Kerrison rongeurs feature a 9mm bite opening and a handle reach 30% shorter than a spurling.

Approximate dimensions:

Jaw Opening (A) - 9 millimeters  
Handle Width (B) - 100 millimeters  
Handle Height (C) - 120 millimeters



### Tarpon™ Style

Tarpon™ from Boss offers the 15mm bite opening of a spurling with the smaller handle size similar to the popular Micro.

Approximate dimensions:

Jaw Opening (A) - 15 millimeters  
Handle Width (B) - 100 millimeters  
Handle Height (C) - 130 millimeters

# APEX KERRISON RONGEURS

PEAK PERFORMANCE AND CLEANABILITY

Introducing APEX Kerrison Rongeurs from BOSS Instruments – one of the most difficult instruments to clean is now one of the easiest!




With a simple *push* of a button, the top rail of the APEX Rongeur lifts up for easy inspection and cleaning.

The top rail locks back in place with a simple squeeze of the handle.

Available in:

- Stainless and Panther coating
- With or without ejector





**Micro Thin Footplate** small handle

**40°↑** 9mm

Length	Style	1.5mm	2mm	3mm	4mm	5mm
* 7" 175mm	Non-Ejector	70-0070	70-0071	70-0072	70-0073	70-0074
	Pin Ejector	--	70-0071E	70-0072E	70-0073E	70-0074E
* 8" 200mm	Non-Ejector	70-0080	70-0082	70-0084	70-0086	70-0088
	Pin Ejector	--	70-0082E	70-0084E	70-0086E	70-0088E

\* Made-to-order. Additional lead time may be required



**APEX Micro Thin Footplate** small handle


**40°↑** 9mm

Length	Style	1mm	2mm	3mm	4mm	5mm
7" 175mm	T-bar Ejector	70-2780TFE	70-2782TFE	70-2784TFE	70-2785TFE	70-2786TFE
8" 200mm	T-bar Ejector	70-2800TFE	70-2802TFE	70-2804TFE	70-2805TFE	70-2806TFE


For Kerrisons with extra performance coatings, please add the following to the end of the item number ( e.g. 70-0200P / 70-0244SC):

**P** = Panther™


**SC** = Super-Coat




Pin ejector



MICRO



T-bar ejector









For Kerrisons with extra performance coatings, please add the following to the end of the item number ( e.g. 70-0200P / 70-0244SC):

**P** = Panther™

**SC** = Super-Coat







## Micro Kerrisons small handle

40°↑ 9mm

Length	Style	1mm 	2mm 	3mm 	4mm 	5mm 	6mm 
7" 175mm	Non-Ejector	70-0160	70-0230	70-0232	70-0234	70-0236	70-0238
	Pin Ejector	--	70-1700	70-1701	70-1702	70-1703	70-1704
8" 200mm	Non-Ejector	70-0162	70-0330	70-0332	70-0334	70-0336	70-0338
	Pin Ejector	--	70-0690	70-0691	70-0692	70-0693	70-0694
9" 225mm	Non-Ejector	70-0164	70-0352	70-0354	70-0356	70-0358	70-0360
	Pin Ejector	--	70-0730	70-0731	70-0732	70-0733	70-0734
10" 250mm	Non-Ejector	70-0166	70-0364	70-0366	70-0368	70-0370	70-0372
	Pin Ejector	--	70-0750	70-0751	70-0752	70-0753	70-0754
11" 275mm	Non-Ejector	70-0168	70-0376	70-0378	70-0380	70-0382	70-0384
	Pin Ejector	--	70-0770	70-0771	70-0772	70-0773	70-0774


## APEX Micro Kerrisons small handle

40°↑ 9mm

Length	Style	1mm 	2mm 	3mm 	4mm 	5mm 	6mm 
7" 175mm	Non-Ejector	70-2780N	70-2782N	70-2784N	70-2785N	70-2786N	70-2787N
	T-bar Ejector	70-2780E	70-2782E	70-2784E	70-2785E	70-2786E	70-2787E
8" 200mm	Non-Ejector	70-2800N	70-2802N	70-2804N	70-2805N	70-2806N	70-2807N
	T-bar Ejector	70-2800E	70-2802E	70-2804E	70-2805E	70-2806E	70-2807E
9" 225mm	Non-Ejector	70-2820N	70-2822N	70-2824N	70-2825N	70-2826N	70-2827N
	T-bar Ejector	70-2820E	70-2822E	70-2824E	70-2825E	70-2826E	70-2827E
*10" 250mm	Non-Ejector	70-2840N	70-2842N	70-2844N	70-2845N	70-2846N	70-2847N
	T-bar Ejector	70-2840E	70-2842E	70-2844E	70-2845E	70-2846E	70-2847E
11" 275mm	Non-Ejector	70-2860N	70-2862N	70-2864N	70-2865N	70-2866N	70-2867N
	T-bar Ejector	70-2860E	70-2862E	70-2864E	70-2865E	70-2866E	70-2867E

\* Made-to-order. Additional lead time may be required

For Kerrisons with extra performance coatings, please add the following to the end of the item number ( e.g. 70-0200P / 70-0244SC ):

 = Panther™

 = Super-Coat



**Ferris-Smith** ergonomic ring handle design

Length	Style	1mm	2mm	3mm	4mm	5mm	6mm
6" (150mm)	Standard	--	70-0810	70-0812	70-0814	70-0816	70-0818
7" (175mm)	Standard	--	70-0834	70-0835	70-0836	70-0837	70-0838
8" (200mm)	Standard	70-0839	70-0840	70-0841	70-0842	70-0843	70-0844
11" (275mm)	Standard	--	70-0910	70-0912	70-0914	70-0916	70-0918

# Boss Instrument Accessories


visit <http://www.bossinst.com> for more information on Boss Instrument accessories.

## Rongeur Cases

Kerrison Rongeur Cases (all hold 8 Rongeurs)  
 00-1220 - 10x12x6" (for up to 9" Rongeurs)  
 00-1230 - 10 x 16 x 6" (for up to 13" Rongeurs)

Foraminotomy Rongeur Case  
 00-1240 - 10 x 12 x 7.5" (extra height for curve)

IVD Rongeur Case (holds up to 16 Rongeurs)  
 00-1250 - 9 x 14 x 5"



00-1220



00-1250

For Kerrisons with extra performance coatings, please add the following to the end of the item number (e.g. 70-0200P / 70-0244SC):







**P** = Panther™

**SC** = Super-Coat



**Foraminotomy** streamlined hingeless design,  
delicate shaft curved up 70°


40°↑ 9mm

Length	Style	1mm 	2mm 	3mm 	4mm 	5mm 	6mm 
8" (200mm)	Standard	--	70-0950	70-0951	70-0952	70-0953	70-0954



**Foraminotomy** hinge design for extra strength,  
shaft available curved up, left and right

40°↑ 9mm

Length	Style	1mm 	2mm 	3mm 	4mm 	5mm 	6mm 
8" (200mm)	Curved Up	--	70-0682	70-0683	70-0684	70-0685	70-0686
	Curved Left	--	--	70-0670	70-0671	70-0672	--
	Curved Right	--	--	70-0674	70-0675	70-0676	--
10" (250mm)	Curved Up	--	70-0682L	70-0683L	70-0684L	--	--

For Kerrisons with extra performance coatings, please add the following to the end of the item number ( e.g. 70-0200P / 70-0244SC):

**P** = Panther™

**SC** = Super-Coat



## IVD Rongeurs


Length	Style	Cushing 1.5mm	Cushing 2mm	Love-Gruenwald 3mm	Spurling 4mm	Spurling 5mm	Cloward 6mm
5" 125mm	Straight	--	70-0980	70-1030	70-1080	--	70-1250
	Up	--	70-0982	70-1032	70-1081	--	70-1251
	Down	--	70-0984	70-1034	70-1082	--	70-1252
6" 150mm	Straight	70-0967	70-0990	70-1040	70-1090	--	70-1430
	Up	70-0968	70-0992	70-1042	70-1091	--	70-1432
	Down	70-0969	70-0994	70-1044	70-1092	--	70-1434
7" 175mm	Straight	70-0970	70-1008	70-1053	70-1100	--	--
	Up	70-0971	70-1009	70-1054	70-1101	--	--
	Down	70-0972	70-1010	70-1055	70-1102	--	--
	Straight w/Teeth	--	70-1120	70-1123	70-1226	--	--
	Up w/Teeth	--	70-1121	70-1124	70-1227	--	--
	Down w/Teeth	--	70-1122	70-1125	70-1228	--	--
9" 225mm	Straight	70-0976	70-1022	70-1070	70-1222	70-1227	--
	Up	70-0977	70-1024	70-1072	70-1224	70-1228	--
	Down	70-0978	70-1026	70-1074	70-1226	70-1229	--
	Straight w/Teeth	70-1021	--	--	--	--	--
10" 250mm	Straight	--	70-1018	70-1066	70-1218	70-1236	--
	Up	--	70-1019	70-1067	70-1219	70-1237	--
	Down	--	70-1020	70-1068	70-1220	70-1238	--
	Straight w/Teeth	--	--	70-1069	70-1221	70-1239	--
11" 275mm	Straight	--	70-1027	70-1076	70-1230	70-1231	--
	Up	--	70-1028	70-1077	70-1232	70-1233	--
	Down	--	70-1029	70-1078	70-1234	70-1235	--
12" 300mm	Straight	--	70-1150	70-1155	70-1160	70-1165	--
	Up	--	70-1151	70-1156	70-1161	70-1166	--
	Down	--	70-1152	70-1157	70-1162	70-1167	--
	Straight w/Teeth	--	--	--	--	70-1239-12	--



## Dragon IVD Rongeurs



Length	Style	2mm x 12mm	3mm x 12mm	4mm x 14mm	5mm x 14mm	6mm x 16mm
5½" 140mm	Straight	70-1281	70-1282	70-1283	--	--
	Up	--	70-1290	70-1291	--	--
	Down	--	70-1296	70-1297	--	--
6¼" 160mm	Straight	70-1285	70-1286	70-1287	70-1288	70-1289
	Up	--	70-1293	70-1294	70-1295	--
	Down	--	70-1298	70-1299	--	--
7¼" 185mm	Straight	70-1256	70-1257	70-1258	70-1259	70-1260
	Up	70-1262	70-1263	70-1264	70-1265	--
	Down	70-1272	70-1273	70-1274	--	--



## Ferris-Smith IVD Rongeurs




Length	Style	2mm	3mm	4mm	5mm	6mm
7" 175mm	Straight	70-1350	70-1351	70-1352	--	70-1354
	Up	70-1360	70-1361	70-1362	--	70-1363
	Down	70-1364	70-1365	70-1366	--	--
10" 250mm	Straight	70-1540	70-1541	70-1542	70-1543	70-1544
	Up	70-1550	70-1551	70-1552	70-1553	70-1554
	Down	70-1555	70-1556	70-1557	70-1558	70-1559
12" 300mm	Straight	--	70-1562	70-1563	70-1564	70-1565
	Up	--	70-1567	70-1568	70-1569	70-1570
	Down	--	70-1572	70-1573	70-1574	70-1575

NOTE: Any Ferris-Smith IVD Rongeur can be ordered With Teeth by adding "WT" to the item number ( e.g. 70-1575WT )

# Modular Kerrison System

All sizes of colour Coded  
Blades fit one handle

Easy push-button blade  
release



The modular Kerrison system saves money in the long term and will always ensure you have sharp Kerrisons for every operation. The interchangeable colour coded blades snap-lock easily into the universal handle, creating a strong instrument. When a blade becomes blunt\*, only a replacement blade\*\* is required, reducing the need to buy a complete new instrument (and reducing spending).

\* when blunt, blades may be resharpened, depending on the condition

\*\* replacement blades are less than 50% the price of a standard Kerrison (April 2015)

8" (20cm) Blade Thin Footplate	8" (20cm) Blade Standard Footplate	9" (23cm) Blade Thin Footplate	9" (23cm) Blade Standard Footplate	Colour Code	Description
PS9000	PS9000	PS9000	PS9000		Modular Handle
PS9011	PS9021	PS9031	PS9041	Red	1mm Blade
PS9012	PS9022	PS9032	PS9042	Yellow	2mm Blade
PS9013	PS9023	PS9033	PS9043	Green	3mm Blade
PS9014	PS9024	PS9034	PS9044	Blue	4mm Blade
PS9015	PS9025	PS9035	PS9045	Black	5mm Blade

- Colour coded blades for easy size identification
- Modular handle fits all blades
- Easy snap-lock /push-button release mechanism
- Blades dismantle for easy cleaning
- Thin & Standard Footplates
- 1 - 5mm bite sizes

We recommend 2 handles per tray and a spare replacement blade for each size (to guarantee a sharp Kerrison for every operation)

