

Osteotomes, Chisels



	Osteotomes	Long Osteotomes	Lateral Bayonet Osteotomes	Thin Shaft Osteotomes	Bayonet Thin Shaft Osteotomes	Sizes
Length	Overall:10.75" (273mm) Handle: 5" (127mm) Working Point:5.75" (146mm)	Overall:15" (380mm) Handle: 5" (127mm) Working Point:10" (254mm)	Overall:18.5" (470mm) Handle: 7" (178mm) Working Point:11.5" (292mm)	Overall:11.5" (292mm) Handle: 4.5" (115MM) Working Point: 7" (178mm)	Overall:11.5" (292mm) Handle: 4.5" (115mm) Working Point: 7" (178mm)	
Straight	PS-O-106	PS-LO-106	PS-LBO-106	PS-TS-106	PS-BTS-106	6mm
Straight	PS-O-109	PS-LO-109	PS-LBO-109	PS-TS-109	PS-BTS-109	9mm
Straight	PS-O-113	PS-LO-113	PS-LBO-113	PS-TS-113	PS-BTS-113	13mm
Straight	PS-O-116	PS-LO-116	PS-LBO-116	PS-TS-116	PS-BTS-116	16mm
Straight	PS-O-119	PS-LO-119	PS-LBO-119	PS-TS-119	PS-BTS-119	19mm
Straight	PS-O-125	PS-LO-125	PS-LBO-125	PS-TS-125	PS-BTS-125	25mm
Curved	PS-O-206	PS-LO-206	PS-LBO-206	PS-TS-206	PS-BTS-206	6mm
Curved	PS-O-209	PS-LO-209	PS-LBO-209	PS-TS-209	PS-BTS-209	9mm
Curved	PS-O-213	PS-LO-213	PS-LBO-213	PS-TS-213	PS-BTS-213	13mm
Curved	PS-O-216	PS-LO-216	PS-LBO-216	PS-TS-216	PS-BTS-216	16mm
Curved	PS-O-219	PS-LO-219	PS-LBO-219	PS-TS-219	PS-BTS-219	19mm
Curved	PS-O-225	PS-LO-225	PS-LBO-225	PS-TS-225	PS-BTS-225	25mm

Material Used: Stainless Steel (AISI-420)

Available as set, or individually, can be provided with customised pattern and sizes.