

Ring Currettes Set



Karlin-Style
Ring Currettes



Straight



6mm



8mm



10mm



Curved



6mm



8mm



10mm

Material Used: Aluminum Handle, Stainless Steel (AISI-420)
Available as set, or individually and Customized Sizes on demand.

	Ring Curette	Karlin-Style Ring Currettes	Sizes
Length	Overall: 12 ³ / ₄ " – 320mm Handle: 7 ³ / ₄ " – 195mm Working Point: 5" – 125mm Handle Diameter: 28mm	Overall: 10" – 260mm Handle: 5" – 130mm Working Point: 5" – 130mm	
Straight	PS-RC-106	PS-KRC-106	6mm
Straight	PS-RC-108	PS-KRC-108	8mm
Straight	PS-RC-110	PS-KRC-110	10mm
Curved	PS-RC-206	PS-KRC-206	6mm
Curved	PS-RC-208	PS-KRC-208	8mm
Curved	PS-RC-210	PS-KRC-210	10mm