


## Spine Colour Currettes



**Material Used:** Handle (Aluminium), Working Points (Titanium / Stainless Steel (AISI-420))

**Length:** 10" (250mm), 4.75" (120mm) working length

Available as set, or individually in titanium or stainless steel

Straight	Forward Angled	Backward Angled	Cup Sizes	Quantity
Overall:10" – 250mm Handle: 5¼" – 135mm Working Point:4¾" – 120mm	Overall:10" – 250mm Handle: 5¼" – 135mm Working Point:4¾" – 120mm	Overall:10" – 250mm Handle: 5¼" – 135mm Working Point:4¾" – 120mm		
PS-S-1405	PS-FA-1405	PS-BA-1405	2	01
PS-S-1406	PS-FA-1406	PS-BA-1406	1	01
PS-S-1411	PS-FA-1411	PS-BA-1411	4-0	01
PS-S-1413	PS-FA-1413	PS-BA-1413	6-0	01