

ULTiMATT™ Colour Blade Cervical Retractor System

PRECISION SURGICAL
Instruments & Retractors of Distinction



The ULTiMATT™ Colour Blade Cervical Retractor System is available in radiolucent Titanium. ULTiMATT™ Black Blade Systems also available

Blade Insertion Clamp





Colour Coded radiolucent Titanium
Blades 22mm and 27mm wide



Radiolucent Lumbar Micro Discectomy Attachments (Extended Set)



ULTiMATT™ Colour Blade Cervical Retractor System



Low Profile, Lightweight
and Strong radiolucent
Titanium Retractor
Frames

Qty	Code	Description
	B23-4000T	ULTiMATT™ Colour Blade TiTANIUM Cervical Retractor System, Colour Blades, Extended
		contains
1	B23-5001T	ULTiMATT™ TiTanium Retractor Frame Hinged Right
1	B23-5002T	ULTiMATT™ TiTanium Retractor Frame Hinged Left
1	B23-5003T	ULTiMATT™ TiTanium Retractor Frame Hinged Offset
1	B23-5006T	ULTiMATT™ TiTanium Long Retractor Rack, 150mm
2	B23-5010	Blade Vertical Insertion Clamps
14	B23-5134T to 70T	Toothed Blades, 22mm wide, 34mm – 70mm, 7 x pairs, Titanium
1	B23-5185T	Toothed Blades, 22mm, 85mm, Titanium
14	B23-5234T to 70T	Blunt Blades, 22mm wide, 34mm – 70mm, 7 x pairs, Titanium
12	B23-5340T to 55T	Two Level Toothed Blades, 27mm wide, 40, 47, 55mm, 4 each size, Titanium
6	B23-5440T to 55T	Two Level Blunt Blades, 27mm wide, 40, 47, 55mm, 3 x pairs, Titanium
1	B23-5030T	Micro Discectomy Attachment, Titanium
1	B23-5021T	Spike Blade, Titanium
1	B23-5023T	T-Blade, Titanium
1	B23-5025T	Interspinous Blade, Titanium
1	B23-5500	Retractor Tray & Lid

Code	Description
B23-4001T	ULTiMATT™ TiTanium Colour Blade, Cervical Retractor System, Standard - same as above (less Long Rack & 4 x Lumbar Discectomy attachments)
B23-5000T	ULTiMATT™ TiTanium Black Blade Titanium, Cervical Retractor System, Black, Extended
B23-5001T	ULTiMATT™ TiTanium Black Blade Cervical Retractor System, Black, Standard same as above (less Long Rack & 4 x Lumbar Discectomy attachments)
B23-4000	ULTiMATT™ Black Blade Cervical Retractor System Black, Extended Stainless Steel
B23-4001	ULTiMATT™ Black Blade Cervical Retractor System, Black, Standard, Stainless Steel same as above (less Long Rack & 4 x Lumbar Discectomy attachments)

To customise your
ULTiMATT™ Colour Blade or
ULTiMATT™ Black Blade
system, please contact us
directly